



# Dispatch® Hospital Cleaner Disinfectant with Bleach

The fastest EPA-registered kill time  
for *Clostridium difficile* spores\*



**Help wipe out MDROs (multidrug-resistant organisms) in just one minute with Dispatch® Hospital Cleaner Disinfectant with Bleach.** Dispatch® spray is a unique, stabilized sodium hypochlorite and detergent solution that kills pathogens fast — with a five-minute contact time for *Clostridium difficile* spores† and one-minute contact time for a broad range of microorganisms, including TB and MDROs such as *Acinetobacter baumannii*, *Enterobacter aerogenes*, VRE and MRSA. (See reverse side for complete efficacy list.)

Dispatch® spray is premixed, ready to use and is similar to the 1:10 bleach solution strength that is recommended by the CDC for effective disinfecting. Dispatch® spray contains an anticorrosive ingredient that makes it less damaging to surfaces or equipment than conventionally prepared bleach solutions.

## Fast-Acting Cleaner Disinfectant

Dispatch® spray is a multipurpose, broad-spectrum cleaner disinfectant for use on hard, nonporous surfaces. It may be used on stainless steel, glazed ceramic tile, hard plastics, nonporous vinyl, painted surfaces, Plexiglas, laminated plastic countertops, fiberglass (sealed) surfaces, glazed porcelain and glass.

- May be used for instrument precleaning and disinfection of noncritical instruments.
- Versatile product may be sprayed directly onto the surface or applied with a clean cloth, mop or sponge saturated with Dispatch® product until the surface is completely wet. A saturated microfiber cloth will deliver optimal surface area coverage.

## Convenient Packaging

- Dispatch® Hospital Cleaner Disinfectant with Bleach is ready to use in convenient spray and refill bottles. No mixing and labeling of bleach solutions is required. A pull-top bottle option for direct application is also available.

\*Based on Federal EPA and State registrations as of September 2010.

†Product with *C. difficile* claim has been registered by the Federal EPA and may not yet be available in all 50 states. Use as directed on precleaned hard nonporous surfaces.

# Dispatch® Hospital Cleaner Disinfectant with Bleach



## Germicidal Efficacy

## KILL TIMES

		KILL TIMES
Bacteria	<i>Acinetobacter baumannii</i>	1 min
	<i>Clostridium difficile</i> spores	5 min
	<i>Enterobacter aerogenes</i>	1 min
	<i>Enterococcus faecium</i> Vancomycin resistant (VRE)	1 min
	<i>Escherichia coli</i> ESBL	1 min
	<i>Klebsiella pneumoniae</i>	1 min
	<i>Pseudomonas aeruginosa</i>	1 min
	<i>Salmonella enterica</i>	1 min
	<i>Staphylococcus aureus</i>	1 min
	<i>Staphylococcus aureus</i> Methicillin resistant (MRSA)	1 min
<i>Streptococcus pyogenes</i>	1 min	
Viruses	Avian Influenza A Virus	1 min
	Canine Parvovirus	1 min
	Feline Panleukopenia Virus	1 min
	Hepatitis A virus (HAV)	1 min
	Hepatitis B virus (HBV)	1 min
	Hepatitis C virus (HCV)	1 min
	Herpes simplex virus (HSV-2)	1 min
	Human Immunodeficiency Virus Type 1 (HIV-1)	1 min
	Influenza A virus	1 min
	Influenza A H1N1	1 min
	Norovirus	1 min
	Poliovirus type 1 (Mahoney)	1 min
	Rhinovirus	1 min
	Rotavirus	1 min
TB	<i>Mycobacterium bovis</i>	30 sec
Fungi	<i>Trichophyton mentagrophytes</i> (Athlete's Foot Fungus)	1 min

## PRODUCT INFORMATION

- Case UPC 68963 15/8 fl. oz. Pump Spray Bottle
- Case UPC 68967 8/22 fl. oz. Trigger Spray Bottle
- Case UPC 68970 6/32 fl. oz. Trigger Spray Bottle
- Case UPC 68832 6/32 fl. oz. Pull-Top Bottle
- Case UPC 68973 6/64 fl. oz. Refill Bottle
- Case UPC 68978 4/128 fl. oz. Refill Bottle
- EPA Reg. No. 56392-7

Dispatch® Hospital Cleaner Disinfectant with Bleach eliminates the labor costs associated with mixing and labeling bleach solutions because it is premixed, ready to use and stable until the expiration date on the package. Unlike bleach and water solutions, Dispatch® spray contains detergent, so it cleans and disinfects in just one step: no precleaning is required.

## Use Sites and Applications

- **Medical:** autoclaves, blood glucose monitors, IV poles, exam tables, bed railings, dialysis machines, respiratory equipment, ultrasound transducers and probes, gurneys, mammography equipment, plastic mattress covers, phlebotomy trays, reception counters, stretchers, wash basins and wheelchairs.
- **Dental:** chairs, counters, operatory hard surfaces and light lens covers.
- **Veterinary:** animal equipment, cages, fountains and clinic surfaces.
- **Laboratory:** countertops, autoclaves and blood glucose monitors.
- **General Use:** doorknobs, hand railings, changing tables, highchairs and playpens.

## Meets Bloodborne Pathogen Standards

OSHA has approved the use of bleach solutions as an effective means of protecting healthcare workers from the risk of contamination by bloodborne pathogens.‡ OSHA requires that conventionally prepared bleach solutions be made up daily (every 24 hours) and left on the surface until the bleach has air dried. Gross contamination (visible soil) must be cleaned up first with soap and water, then followed by use of the bleach solution for disinfecting.

‡Source: OSHA Instruction CPL 02-02-069, November 27, 2001.

Use as directed on hard, nonporous surfaces.

## Expiration Date

The product expiration date should be visible on the package. If not, contact us at [healthcare@clorox.com](mailto:healthcare@clorox.com).

## Safety Information

**Material Safety Data Sheets:** Complete MSDS information can be found on [www.cloroxprofessional.com](http://www.cloroxprofessional.com).

**Precautions:** Dispatch® spray should be used in a ventilated area — general room ventilation is adequate for most users.

For more information, contact your Clorox sales representative or call **800-234-7700**.

email: [healthcare@clorox.com](mailto:healthcare@clorox.com)

visit us: [cloroxprofessional.com](http://cloroxprofessional.com)

