



HEALTHCARE SOLUTIONS

CLOROX PROFESSIONAL PRODUCTS

CLOROX

Clorox® Germicidal Wipes

- Ready-to-use, hospital-grade disinfectant
- Kill claims for 36 different microorganisms, including *Clostridium difficile* spores, MRSA, HIV-1, *Pseudomonas aeruginosa*, *Staphylococcus aureus*, *Streptococcus pyogenes*, *Streptococcus pneumoniae* and *Acinetobacter baumannii*
- EPA-registered to kill *Clostridium difficile* spores in 5 minutes*
- Contact time of 30 seconds for most bacteria (see product label)
- Contains the power of Clorox® Bleach
- 0.55% sodium hypochlorite/5500 ppm
- Convenient wipe application
- EPA Reg. No. 67619-12



The fastest EPA-registered kill time for *Clostridium difficile* spores*

Clorox® Germicidal Wipes, Terminal
12" x 12" wipe size

Case UPC 30358 2/110 ct.
Case UPC 30359 2/110 ct. Refill

Clorox® Germicidal Wipes, Clinical
6" x 5" wipe size

Case UPC 30577 6/150 ct.

Clorox® Germicidal Wipes, Multi-Purpose
6.75" x 9" wipe size

Case UPC 35309 6/70 ct.

Dispatch® Hospital Cleaner Disinfectant with Bleach

- Premixed, ready-to-use 1:10 bleach solution at a concentration that is recommended by CDC guidelines for effective disinfecting
- EPA-registered to kill *Clostridium difficile* spores in 5 minutes*
- 1-minute contact time for bacteria, viruses and fungi, and 30 seconds for TB
- Cleans and disinfects in just one step; no precleaning required
- Available in ready-to-use spray for broad coverage, pull-top for direction application and convenient refill bottles
- EPA Reg. No. 56392-7



The fastest EPA-registered kill time for *Clostridium difficile* spores*

Dispatch® Hospital Cleaner Disinfectant with Bleach, Spray

Case UPC 68970 6/32 fl. oz.

Dispatch® Hospital Cleaner Disinfectant with Bleach, Pull-top

Case UPC 68832 6/32 fl. oz.

Dispatch® Hospital Cleaner Disinfectant with Bleach, Refill

Case UPC 68978 4/128 fl. oz.

NEW!

Clorox Commercial Solutions® Clorox® Bleach Cream Cleanser

Case UPC 30613 8/32 fl. oz.

- Tough cleaning power leaves surfaces sanitary
- Commercial-grade cleaning
- Removes stubborn stains like coffee, juice, soap scum and blood
- Easy rinse; clean, bright shine; gentle on surfaces
- Color shows coverage for easy, complete cleaning



Dispatch® Hospital Cleaner Disinfectant with Bleach

Case UPC 68967 8/22 fl. oz. Spray

Case UPC 68973 6/64 fl. oz. Refill

EPA Reg. No. 56392-7

- Premixed, ready-to-use solution similar to 1:10 bleach solution strength recommended by CDC guidelines for effective disinfecting
- EPA-registered to kill *Clostridium difficile* spores in 5 minutes*
- 1-minute contact time for bacteria, viruses and fungi, and 30 seconds for TB
- Cleans and disinfects in just one step; no precleaning required
- Also available in 22 oz. and 64 oz. refill

The fastest
EPA-registered
kill time for
Clostridium difficile
spores*



Dispatch® Hospital Cleaner Disinfectant Towels with Bleach

Case UPC 69101 6/50 ct. (7" x 8" wipe size), Individually Wrapped

Case UPC 69240 24/40 ct. (8" x 10" wipe size), Patient Transport Pack

Case UPC 69260 12/60 ct. (9" x 10" wipe size), Soft Pack

Case UPC 69150 8/150 ct. (6.75" x 8" wipe size), Canister

EPA Reg. No. 56392-8

- Premoistened towels with a 1:10 bleach solution strength recommended by CDC guidelines for effective disinfecting
- EPA-registered to kill *Clostridium difficile* spores in 5 minutes*
- 1-minute contact time for a broad range of organisms, including MDROs, and 2 minutes for TB
- Cleans and disinfects in just one step; no precleaning required
- Available in transport pack for mobile equipment, soft pack for carts and individual for on-the-go needs

The fastest
EPA-registered
kill time for
Clostridium difficile
spores*



Clorox Commercial Solutions® Ultra Clorox® Germicidal Bleach

Case UPC 02490 6/96 fl. oz.

Case UPC 02489 3/182 fl. oz.

EPA Reg. No. 67619-8

- The first germicidal bleach registered by the EPA to kill *C. difficile* spores†
- Kill claims for 38 different microorganisms, including *Clostridium difficile* spores, TB, HIV-1, Hepatitis A, B and C viruses, MRSA, VRE, *Streptococcus pneumoniae*, Human Coronavirus, Norovirus (as feline calicivirus), *Acinetobacter baumannii* and Avian influenza A
- 6.15% sodium hypochlorite/61500 ppm
- Helps with compliance with OSHA bloodborne pathogen guidelines



*Based on federal EPA registrations as of July 2011 (5-minute contact time). *C. difficile* spore claim has been registered by the federal EPA and may not yet be available in all 50 states. Check with your sales representative for updates in your state. Use as directed on pre-cleaned hard, nonporous surfaces.

†*C. difficile* claim federal EPA registered in December 2008.

Non-Bleach-Based Products

Citrex® Hospital Spray Disinfectant

Case UPC 49200 12/14 oz.
EPA Reg. No. 56392-4

- 1-minute kill time for TB, HIV-1, Hepatitis B, Hepatitis C
- 2-minute kill time for *Staphylococcus aureus*, *Pseudomonas aeruginosa*, *Salmonella choleraesuis*, *Streptococcus pyogenes*
- Convenient, easy to use and economical
- Designed for use in all types of medical facilities



Citrace® Hospital Germicide

Case UPC 49100 12/14 oz.
EPA Reg. No. 56392-2

- Eliminates odors
- Kills 99.9% of bacteria on hard surfaces in 30 seconds
- Disinfects hard surfaces in 5 minutes or less
- Inhibits the growth of mold and mildew on fabric



Precise® Hospital Foam Cleaner Disinfectant

Case UPC 40513 12/20 oz.
EPA Reg. No. 56392-1

- Cleans and disinfects in one step; no surface precleaning required
- 1-minute contact time for Hepatitis B, Hepatitis C and HIV-1
- 10-minute contact time for hospital bacteria, TB and fungi
- Use on wide range of surfaces and medical/dental equipment



Precise® QTB Hospital Cleaner Disinfectant

Case UPC 23032 9/32 fl. oz., Spray
Case UPC 23128 4/128 fl. oz., Refill
EPA Reg. No. 1839-83-56392

- One-step, ready-to-use cleaner disinfectant
- Low-toxicity quaternary/water formulation can be used on washable surfaces such as metal, stainless steel, glazed porcelain, glazed ceramic tile and plastic
- Disinfects food service establishment food-contact surfaces
- **In 30 Seconds:** kills Feline Calicivirus, Norovirus, Rabies Virus
- **In 1 Minute:** kills Human Immunodeficiency Virus Type 1 (HIV-1)
- **In 2 Minutes:** kills Avian Influenza A Strain H3N2, Avian Influenza A Strain H9N2, Human Coronavirus, SARS-associated Coronavirus
- **In 3 Minutes:** kills hospital bacteria, other bacteria and viruses
- **In 5 Minutes:** kills TB



Clorox Commercial Solutions® Clorox® Disinfecting Spray

Case UPC 38504 12/19 oz.
EPA Reg. No. 67619-3

- Sanitizes, disinfects and deodorizes
- Excellent for use in nursing homes, hospitals, day care centers, schools and hotels
- Effective quat-based and ethanol-based disinfectant against E. coli, TB, HIV-1, Hepatitis A virus and MRSA
- Fresh scent
- Bleach-free



Clorox Commercial Solutions® Clorox® Hand Sanitizer

Spray Dispenser, Case UPC 01752 6/1 L
Spray Refill, Case UPC 01753 6/1 L
Spray, Case UPC 02176 12/500 mL
Spray, Case UPC 02174 24/2 oz.

- Convenient spray form to sanitize hands
- 71% ethyl alcohol
- Effective against 77 microorganisms
- Kills 99.999% of TFM required germs
- Moisturizes hands without sticky residue
- Ideal for use when soap and water are not available



For more information, contact your Clorox sales representative or call **(800) 234-7700**.
email: healthcare@clorox.com
visit us: cloroxprofessional.com

© 2011 Clorox Professional Products Company,
1221 Broadway, Oakland, CA 94612